



# ARIAS LABS SYSTEM

## FURNITURE & EQUIPMENT

SAFETY  
QUALITY  
ADAPTABILITY



THINK  
LABS  
THINK  
ARIAS



Calle N3, Parcela C  
Polígono industrial el Bobalar,  
46130 Massamagrell, Valencia  
Tel: +34 960 013 119 | +34 634 652 807  
[www.ariaslabssystem.com](http://www.ariaslabssystem.com)

# WE ARE

---

- Laboratory Furniture Designers & Manufacturers
- Hospital Furniture Solution Providers

# ARIAS

**Accessibility**

**Reliability**

**Innovation**

**Adaptability**

**Safety**







## ABOUT US



**ARIAS LABS SYSTEM** is a Spanish laboratory furniture, fume hood and laboratory equipment manufacturer with more than 20 years of experience in this field.

Arias is a complete solution provider for your laboratory needs. Our services include designing, planning, engineering, commissioning & consulting for all types of laboratories.

We provide our clients with the best solutions required for their laboratory function, having in mind the best proposition in terms of safety, quality and delivery time.

For different target groups & applications, our ergonomic design enables researchers to work comfortably and safely in any research environment.



## OUR VISION



Getting a safe laboratory with innovative design should be accessible to any company or organization, in terms of budget, product adaptability, right design & delivery time. That is where we focus while designing, engineering & manufacturing new products. Arias aims to be one of the best innovative European laboratory manufacturers operating worldwide.

# INDEX

- **Material Specifications and Manufacturing Process ..... 07**
  
- **Our Products ..... 11**
  - Fume Hoods ..... 12
  - Biological Safety Cabinets ..... 24
  - Laboratory Benches ..... 32
  - Worktops ..... 41
  - Sink Units ..... 43
  - Storage Units ..... 46
  - Service Galleries ..... 59
  - Wall Partitions & Service Column ..... 67
  - Accessories ..... 70
  
- **Product Reference List ..... 75**
  
- **Certifications ..... 83**



# MATERIAL SPECIFICATIONS AND MANUFACTURING PROCESS








## MATERIAL SPECIFICATIONS AND MANUFACTURING PROCESS

### STEEL SHEET

Our bench frames, steel cabinets & cupboards are manufactured by electro-zinc cold rolled steel epoxy polyester coated resistant to corrosion. Before painting, the steel sheets go through: electrolytic zinc settlement, degreasing & surface cleaning, pre-treatment with phosphate and chromated treatment.



Particle Board

MDF Melamine



### MELAMINE SHEET

For some specific projects, a line of cabinets & cupboards is manufactured from wood fibres and agglutinated by synthetic resins with pressure and temperature. The MDF melamine is moisture resistant, with low formaldehyde content, compact boards and lightweight.



## PAINTING PROCESS

Our products are coated by a chemical resistant epoxy powder paint. The painting is specially formulated to satisfy the laboratory needs in order to obtain a very high resistance to UV radiations, chemicals and to humidity. Its excellent gloss retention and resistance to colour change, combined with other outstanding properties are what makes Arias different.

### **Our painting process follows these meticulous steps:**

- 1- Surface polishing & dedusting if needed.
- 2- Surface cleaning: Acid washing & phosphating or nano zirconia treatment.
- 3- Washing, drying & surface examination.
- 4- Acid resistant epoxy powder automatic painting.
- 5- Drying at 200°C in a special oven.
- 6- Cooling & packing.





A long, white laboratory bench with multiple drawers and cabinets, set against a white wall with a grid ceiling. The bench is equipped with several drawers and cabinets, some of which are on wheels. The text "OUR PRODUCTS" is overlaid in the center of the image.

## OUR PRODUCTS





## FUME HOODS





Fume Hoods are ventilated enclosures that enable laboratory personnel to safely work with materials that produce hazardous fumes and volatile vapors.

We at Arias Labs continue to innovate the standards of our fume hoods to support the evolving needs of modern laboratory safety requirements & end user new application requirements.

Our fume hoods are suitable for general and special laboratory applications, dealing with various common acid and solvents as well as enduring the harshest laboratory environments. Our fume hoods are certified as per standards EN14175 & ASHRAE 110-2016.

## Fume Hood Types:

- Basic Fume Hood – Manual Control Panel
- Standard Fume Hood – LCD Control Panel
- Prime Fume Hood – LCD Control Panel Energy Saving Inverter (Automatic Flow)
- Walk-in Fume Hood
- Polypropylene Fume Hood – Perchloric HCl Acid FH





### Fume Hood Basic Type - Manual Control Panel

This fume hood is suitable for general laboratory applications, dealing with various common acid and solvents as well as enduring the harshest laboratory environments.

Fume hood basic type could be used as a ductless fume hood as well. This type of fume hood utilizes carbon or HEPA filtration to protect laboratory personnel from toxic chemical fumes, odors, and particles. Ductless Fume Hood guarantees a hassle-free installation and no lab work interruption.





## FUME HOODS

Model	AR01-120	AR01-150	AR01-180	AR01-210
<b>External dimensions (mm) (W x D x H)</b>	1200x750x2200	1500x750x2200	1800x750x2200	2100x750x2200
<b>Upper Body</b>	Steel epoxy powder coated for chemical resistance			
<b>Support Stand</b>	Integrated (full height body)			
<b>Color</b>	Light Grey & Dark Grey			
<b>Worktop</b>	16 mm Anti-chemical compact laminate countertop			
<b>Under FH cabinet</b>	Vented Cabinet steel epoxy powder coated			
<b>Operation Controller Panel</b>	<ul style="list-style-type: none"> <li>• Main fuse breaker</li> <li>• Manual Fan control</li> <li>• Manual Speed control</li> <li>• Manual Lamp switch</li> </ul>			
<b>Sash type</b>	Vertical opening			
<b>Sash material</b>	Tempered safety glass 6mm with aluminum frame			
<b>Sash Safety door stop</b>	Yes			
<b>Lamp</b>	LED 1 unit			
<b>Cup sink</b>	PP 1 unit			
<b>Water tap</b>	CW 1 unit			
<b>Air valve</b>	Option			
<b>Vacuum valve</b>	Option			
<b>Sockets</b>	British type 2 units with cover (international type available)			
<b>Electrical Break</b>	16A Insurance			
<b>Duct collar (Ø)</b>	PP 250mm (Ø)			
<b>Duct</b>	PVC flexible duct (3 m each fume hood)			
<b>Blower Type: Polypropylene built in blower</b>	1040 m³/h	1400 m³/h	2350 m³/h	3400 m³/h
<b>Option</b>	Worktop: Ceramic - Color: Off white & Light Grey Controller: Digital touch panel - Sash: Horizontal opening Air, Vacuum Valve			



## Fume Hood Standard Type - LCD Control Panel

This fume hood is suitable for general laboratory applications, dealing with various common acid and solvents as well as enduring the harshest laboratory environments.

This fume hood comes with LCD Display control panel for AIRFLOW controlling and monitoring, in addition to visual & audible alarms, air flow digital screen, and touch sensitive buttons for all the operations.





## FUME HOODS

Model	AR02-120	AR02-150	AR02-180	AR02-210
<b>External dimensions (mm) (W x D x H)</b>	1200x900x2340	1500x900x2340	1800x900x2340	2100x900x2340
<b>Upper Body</b>	Steel epoxy powder coated for chemical resistance			
<b>Support Stand</b>	Pipe with horizontal supports epoxy powder coated			
<b>Inner services connection</b>	From the back side of the Fume hood			
<b>Frontal window</b>	700mm			
<b>Color</b>	Light Grey & Dark Grey			
<b>Worktop</b>	16 mm Anti-chemical compact laminate countertop			
<b>Under FH cabinet</b>	Vented Cabinet steel epoxy powder coated			
<b>Operation Controller Panel</b>	LCD Display control panel for AIRFLOW controlling and monitoring <ul style="list-style-type: none"> <li>• Visual &amp; audible alarms</li> <li>• Air flow digital screen</li> <li>• Touch sensitive buttons for all the operations</li> </ul>			
<b>Sash type</b>	Vertical opening			
<b>Sash material</b>	Tempered safety glass 6mm with aluminum frame			
<b>Sash Safety door stop</b>	Steel epoxy powder coated			
<b>Lamp</b>	LED 1 unit			
<b>Cup sink</b>	PP 1 unit			
<b>Water tap</b>	CW 1 unit			
<b>Air valve</b>	CA 1 unit			
<b>Vacuum valve</b>	Option			
<b>Sockets</b>	British type 2 units with cover (international type available)			
<b>Electrical Break</b>	16A Insurance			
<b>Duct collar (Ø)</b>	PP 250mm (Ø)			
<b>Duct</b>	PVC flexible duct (3 m each fume hood)			
<b>Blower Type: Polypropylene built in blower</b>	1040 m³/h	1400 m³/h	2350 m³/h	3400 m³/h
<b>Option</b>	Worktop: Ceramic – Cup sink: Ceramic Color: Off white & Light Grey – Sash: Horizontal opening Fittings: Vacuum Valve			





### **Fume Hood Prime Type – LCD Control Panel Energy Saving Inverter (Auto- matic Flow)**

This fume hood is suitable for general & special laboratory applications, dealing with various common acid and solvents as well as enduring the harshest laboratory environments.

This fume hood comes with VAC (variable air control) system which enables the control the speed of the exhaust air to save energy & fresh air. Also, it comes with LCD Display control panel for AIRFLOW controlling and monitoring, in addition to visual & audible alarms, air flow digital screen, and touch sensitive buttons for all the operations.





## FUME HOODS

Model	AR03-120	AR03-150	AR03-180	AR03-210
<b>External dimensions (mm) (W x D x H)</b>	1200x817x2380	1500x817x2380	1800x817x2380	2100x817x2380
<b>Upper Body</b>	Steel epoxy powder coated for chemical resistance			
<b>Support Stand</b>	Pipe with horizontal supports epoxy powder coated			
<b>Inner services connection</b>	From the lateral sides of the fume hood			
<b>Frontal window</b>	1100mm			
<b>Color</b>	Light Grey & Dark Grey			
<b>Worktop</b>	Ceramic worktop with 28mm thick edge marine			
<b>Under FH cabinet</b>	Vented Cabinet steel epoxy powder coated			
<b>Operation Controller Panel</b>	LCD Display control panel for AIRFLOW controlling and monitoring <ul style="list-style-type: none"> <li>• Visual &amp; audible alarms</li> <li>• Air flow digital screen</li> <li>• Touch sensitive buttons for all the operations</li> </ul>			
<b>Sash type</b>	Vertical opening			
<b>Sash material</b>	Tempered safety glass 6mm with aluminum frame			
<b>Sash Safety door stop</b>	Steel epoxy powder coated			
<b>Lamp</b>	LED 1 unit			
<b>Cup sink</b>	PP 1 unit			
<b>Water tap</b>	CW 1 unit			
<b>Air valve</b>	CA 1 unit			
<b>Vacuum valve</b>	Option			
<b>Sockets</b>	British type 2 units with cover (international type available)			
<b>Electrical Break</b>	16A Insurance			
<b>Duct collar (Ø)</b>	PP 250mm (Ø)			
<b>Duct</b>	PVC flexible duct (3 m each fume hood)			
<b>Blower Type: Polypropylene built in blower</b>	1040 m <sup>3</sup> /h	1400 m <sup>3</sup> /h	2350 m <sup>3</sup> /h	3400 m <sup>3</sup> /h
<b>Option</b>	Worktop: Ceramic – Cup sink: Ceramic Color: Off white & Light Grey – Sash: Horizontal opening Fittings: Vacuum Valve			



## Walk-in Fume Hood

Walk-In Fume Hoods are "oversized" enclosures that are built to isolate the exterior of the hood from fumes, dust, odors, and other harmful particulate. Primarily, these hoods isolate the environment from particle contamination and allow applications to be performed in a sperate working area.





## FUME HOODS

Model	ARO4-120	ARO4-150	ARO4-180	ARO4-210
<b>External dimensions (mm) (W x D x H)</b>	1200x900x2340	1500x900x2340	1800x900x2340	2100x900x2340
<b>Upper Body</b>	Steel epoxy powder coated for chemical resistance			
<b>Support Stand</b>	Pipe with horizontal supports epoxy powder coated			
<b>Inner walls</b>	One side to allow the connections to be in the back side of the Fume hood			
<b>Frontal window</b>	700mm opening space for the operations			
<b>Color</b>	Light Grey & Dark Grey			
<b>Worktop</b>	N/A			
<b>Under FH cabinet</b>	N/A			
<b>Operation Controller Panel</b>	LCD Display control panel for AIRFLOW controlling and monitoring <ul style="list-style-type: none"> <li>• Visual &amp; audible alarms</li> <li>• Air flow digital screen</li> <li>• Touch sensitive buttons for all the operations</li> </ul>			
<b>Sash type</b>	Vertical opening			
<b>Sash material</b>	Double Tempered safety glass 6mm with aluminum frame			
<b>Sash Safety door stop</b>	Steel epoxy powder coated			
<b>Lamp</b>	LED 1 unit			
<b>Cup sink</b>	PP 1 unit			
<b>Water tap</b>	CW 1 unit			
<b>Air valve</b>	CA 1 unit			
<b>Vacuum valve</b>	Option			
<b>Sockets</b>	British type 2 units with cover (international type available)			
<b>Electrical Break</b>	16A Insurance			
<b>Duct collar (Ø)</b>	PP 250mm (Ø)			
<b>Duct</b>	PVC flexible duct (3 m each fume hood)			
<b>Blower Type: Polypropylene built in blower</b>	1040 m³/h	1400 m³/h	2350 m³/h	3400 m³/h
<b>Option</b>	Cup sink: Ceramic Color: Off white & Light Grey Fittings: Vacuum Valve			



### **Polypropylene Fume Hood – Perchloric HCl Acid FH**

Polypropylene Fume Hoods are made of high-quality non-corrosive polypropylene with excellent chemical resistance. Polypropylene increases the product's tensile strength and improves its thermal characteristics. The Polypropylene Fume Hoods protect laboratory staff from noxious fumes released by acids, dangerous gas and organic solutions – materials and acids which regular steel hoods may not withstand.







Model	AR05-120	AR05-150	AR05-180	AR05-210
<b>External dimensions (mm) (W x D x H)</b>	1200x817x2380	1500x817x2380	1800x817x2380	2100x817x2380
<b>Upper Body</b>	One-piece welded Polypropylene structure			
<b>Support Stand</b>	Integrated (full height body)			
<b>Color</b>	Beige			
<b>Worktop</b>	Built-in polypropylene worktop			
<b>Under FH cabinet</b>	Integrated Vented Cabinet full polypropylene			
<b>Operation Controller Panel</b>	LCD Display control panel for AIRFLOW controlling and monitoring <ul style="list-style-type: none"> <li>• Visual &amp; audible alarms</li> <li>• Air flow digital screen</li> <li>• Touch sensitive buttons for all the operations</li> </ul>			
<b>Sash type</b>	Vertical opening			
<b>Sash material</b>	Tempered safety glass 6mm			
<b>Lamp</b>	LED 1 unit			
<b>Cup sink</b>	PP 1 unit			
<b>Water tap</b>	CW 1 unit			
<b>Air valve</b>	CA 1 unit			
<b>Vacuum valve</b>	Option			
<b>Sockets</b>	British type 2 units with cover (international type available)			
<b>Electrical Break</b>	16A Insurance			
<b>Duct collar (Ø)</b>	PP 250mm (Ø)			
<b>Duct</b>	PVC flexible duct (3 m each fume hood)			
<b>Blower</b>	Option			
<b>Option</b>	Worktop: Ceramic – Cup sink: Ceramic Sash: Horizontal opening – Fittings: Vacuum Valve Blower: Centrifugal FAN			



## BIOLOGICAL SAFETY CABINETS



## BIOLOGICAL SAFETY CABINETS

A Biological Safety Cabinet is mainly used for handling pathogenic biological samples or for applications that require a sterile work zone. It creates inflow and downflow of air that ensures the user, product and environment protection from the biohazard factor and prevents cross-contamination of various samples used in experiments.

### Biosafety Cabinets Types:

- Class I: Laminar Flow
- Class II: A2 Biological Safety Cabinet
- Class II: B2 Biological Safety Cabinet





## Class I: Laminar Flow

A laminar flow cabinet is defined as enclosed workbench which is used to create a contamination free work environment through installed HEPA filters that capture all the particles entering the cabinet.

A laminar flow cabinet is used for work with substances which are not hazardous for the personnel health for personnel and environmental protection.

### Laminar Flow Types:

- Horizontal Laminar Flow
- Vertical Laminar Flow
- PCR Workstation





## BIOLOGICAL SAFETY CABINETS

Model	AR-LAM-90	AR-LAM-120	AR-LAM-150	AR-LAM-180
<b>External Dimensions (mm) (W x D x H)</b>	1065x900x1550	1360x900x1550	1660x900x1550	1960x900x1550
<b>Internal Dimensions (mm) (W x D x H)</b>	900x630x650	1200x630x650	1500x630x650	1800x630x650
<b>Stand Dimensions (mm) (W x D x H)</b>	1060x870x850	1360x870x850	1660x870x850	1960x870x850
<b>Exterior Body</b>	Steel epoxy powder coated for chemical resistance			
<b>Inner walls</b>	Stainless steel (AISI 304)			
<b>Worktop</b>	Worktop floor height 850mm made of Stainless steel (AISI 316)			
<b>Operation Controller Panel</b>	LCD Display control panel for AIRFLOW controlling and monitoring <ul style="list-style-type: none"> <li>• Visual &amp; audible alarms</li> <li>• Air flow digital screen</li> <li>• Touch sensitive buttons for all the operations</li> </ul>			
<b>Average Inflow Velocity m/s</b>	0.45			
<b>Average Downflow Velocity m/s</b>	0,25 / 0,50			
<b>Filters (EN1822) Main Filter</b>	H14 HEPA, 0.3 µm particulate 99.995%			
<b>Filters (EN1822) Exhaust Filter</b>	Not used in Laminar flow Cabinet			
<b>Filters (EN1822) Front Filter</b>	G4			
<b>Cleanliness level (US FED 209E)</b>	ISO 5 - Class 100			
<b>Exhaust Fan Motor Capacity</b>	Suitable for 10 m Horizontal and 3 m Vertical Long Air Channel			
<b>Front Sash / Window Standard Opening</b>	200mm ± 10mm			
<b>Front Sash/Window Maximum Opening</b>	530mm ± 10mm			
<b>Noise Level (Normal Mode)</b>	< 65 dB(A)			
<b>Noise Level (Economy Mode)</b>	< 55 dB(A)			
<b>Illumination</b>	850 – 1150 Lux			
<b>Front Sash/Window Thickness</b>	< 6 mm			
<b>Power Consumption (Fan Filter Unit)</b>	105 W		200 W	
<b>Power Consumption (Lighting)</b>	50 W			
<b>Power Consumption (UV Lamp)</b>	30 W			
<b>Power supply (Voltage) Socket</b>	220-240 V / 16A			
<b>Supply Voltage / Frequency / Power / Current</b>	230 V / 50 Hz / 2500 W / 16A			
<b>Cabin Weight</b>	205 kg	230 kg	270 kg	285 kg
<b>Option</b>	Fittings: Vacuum Valve			



### Class II: A2 Biological Safety Cabinet

Class II: A2 Biological Safety Cabinet ensures the user, product and environment protection from the biohazard factor and prevents cross-contamination of various samples used in experiments. They recirculate 70% of the HEPA-filtered laminar air while 30% passes through an exhaust HEPA filter for discharge to the lab or optional canopy.







## BIOLOGICAL SAFETY CABINETS

Model	ARS-2A2-90	ARS-2A2-120	ARS-2A2-150	ARS-2A2-180
<b>External Dimensions (mm) (W x D x H)</b>	1065x900x1450	1360x900x1450	1660x900x1450	1960x900x1450
<b>Internal Dimensions (mm) (W x D x H)</b>	900x630x650	1200x630x650	1500x630x650	1800x630x650
<b>Stand Dimensions (mm) (W x D x H)</b>	1060x870x750	1360x870x750	1660x870x750	1960x870x750
<b>Exterior Body</b>	Steel epoxy powder coated for chemical resistance			
<b>Inner walls</b>	Stainless steel (AISI 304)			
<b>Worktop</b>	Worktop floor height 850mm made of Stainless steel (AISI 316)			
<b>Operation Controller Panel</b>	LCD Display control panel for AIRFLOW controlling and monitoring <ul style="list-style-type: none"> <li>• Visual &amp; audible alarms</li> <li>• Air flow digital screen</li> <li>• Touch sensitive buttons for all the operations</li> </ul>			
<b>Average Inflow Velocity m/s</b>	0.45			
<b>Average Downflow Velocity m/s</b>	0,25 / 0,50			
<b>Filters (EN1822) Main Filter</b>	H14 HEPA, 0.3 µm particulate 99.995%			
<b>Filters (EN1822) Exhaust Filter</b>	H14 HEPA, 0.3 µm particulate 99.995%			
<b>Filters (EN1822) Front Filter</b>	Not used in Class II A II			
<b>Cleanliness level (US FED 209E)</b>	ISO 5 - Class 100			
<b>Exhaust Fan Motor Capacity</b>	Air flow Qv 2400 m³/h			
<b>Front Sash / Window Standard Opening</b>	200mm ± 10mm			
<b>Front Sash/Window Maximum Opening</b>	530mm ± 10mm			
<b>Noise Level (Normal Mode)</b>	< 65 dB(A)			
<b>Noise Level (Economy Mode)</b>	< 55 dB(A)			
<b>Illumination</b>	850 – 1150 Lux			
<b>Front Sash/Window Thickness</b>	< 6 mm			
<b>Power Consumption (Fan Filter Unit)</b>	105 W		200 W	
<b>Power Consumption (Lighting)</b>	50 W			
<b>Power Consumption (UV Lamp)</b>	30 W			
<b>Power supply (Voltage) Socket</b>	220-240 V / 16A			
<b>Supply Voltage / Frequency / Power / Current</b>	230 V / 50 Hz / 2500 W / 16A			
<b>Cabin Weight</b>	205 kg	230 kg	270 kg	285 kg
<b>Option</b>	Fittings: Vacuum Valve			



### Class II: B2 Biological Safety Cabinet

Class II: B2 Biological Safety Cabinet ensures the user, product and environment protection from the biohazard factor and prevents cross-contamination of various samples used in experiments. The type B2 total exhaust cabinets are widely used in toxicology laboratories and similar applications where chemical effluent is present and clean air is essential.

The air goes through a HEPA filter and is 100% exhausted through a dedicated duct with an exhaust fan motor. The air is then released into the atmosphere.





## BIOLOGICAL SAFETY CABINETS

Model	ARS-2B2-90	ARS-2B2-120	ARS-2B2-150	ARS-2B2-180
External Dimensions (mm) (W x D x H)	1065x900x1450	1360x900x1450	1660x900x1450	1960x900x1450
Internal Dimensions (mm) (W x D x H)	900x630x650	1200x630x650	1500x630x650	1800x630x650
Stand Dimensions (mm) (W x D x H)	1060x870x750	1360x870x750	1660x870x750	1960x870x750
Exterior Body	Steel epoxy powder coated for chemical resistance			
Inner walls	Stainless steel (AISI 304)			
Worktop	Worktop floor height 850mm made of Stainless steel (AISI 316)			
Operation Controller Panel	LCD Display control panel for AIRFLOW controlling and monitoring <ul style="list-style-type: none"><li>• Visual &amp; audible alarms</li><li>• Air flow digital screen</li><li>• Touch sensitive buttons for all the operations</li></ul>			
Average Inflow Velocity m/s	0.45			
Average Downflow Velocity m/s	0,25 / 0,50			
Filters (EN1822) Main Filter	H14 HEPA, 0.3 µm particulate 99.995%			
Filters (EN1822) Exhaust Filter	H14 HEPA, 0.3 µm particulate 99.995%			
Filters (EN1822) Front Filter	G4			
Cleanliness level (US FED 209E)	ISO 5 - Class 100			
Exhaust Fan Motor Capacity	Air flow Qv 2400 m³/h			
Front Sash / Window Standard Opening	200mm ± 10mm			
Front Sash / Window Maximum Opening	530mm ± 10mm			
Noise Level (Normal Mode)	< 65 dB(A)			
Noise Level (Economy Mode)	< 55 dB(A)			
Illumination	850 – 1150 Lux			
Front Sash / Window Thickness	6 mm			
Power Consumption (Fan Filter Unit)	105 W		200 W	
Power Consumption (Lighting)	50 W			
Power Consumption (UV Lamp)	30 W			
Power supply (Voltage) Socket	220-240 V / 16A			
Supply Voltage / Frequency / Power / Current	230 V / 50 Hz / 2500 W / 16A			
Weight	205 kg	230 kg	270 kg	285 kg
Option	Fittings: Vacuum Valve, gases			



LABORATORY BENCHES



**A wide range of modular configurations are available for any type of laboratory**

### TABLE CONFIGURATION

Arias Labs System provides modular solutions for any special requirement and for all challenging space constraints with its versatile design. All our benches can be selected with various worktop materials for a variety of application possibilities in the laboratory. They are designed in a way that ensures flexibility and high levels of stability in order to provide you with the perfect platform for your research & equipment.

**Arias offers 3 types of bench systems:**

- C-frame system
- H-frame system
- On-plinth system

All our benches' frames are manufactured of epoxy polyester coated product that are made of zinc coated cold rolled steel resistant to chemical reagents.



Available Modular Configurations:

### C-FRAME SYSTEM

The C-frame system was designed to provide a flexible modular solution with under bench registrable cover panels for easy maintenance & disinfection. It can incorporate different types of under bench cabinets: suspended, sliding or on-wheels depending on the client's needs. Our c-frames are extremely stable and can load up to 200 kg/m<sup>2</sup>.

C-frames are designed for different positions and styles:

- Wall bench
- Island bench
- Movable bench







## LABORATORY BENCHES


### Wall Laboratory Bench:

These frames are provided with a leveling system in contact with the floor as well as in contact with the working surfaces.

On the rear, it provides a structural support for the fixation of the service system, separated from the frame of the workbench, creating a space (service gallery) for the guiding of the mechanical, electrical, plumbing, and connections, easily accessible via removable covers.

Additional shelves could be added upon the request of the clients. These storage shelves are anchored to the vertical profiles and easily height adjustable. They are protected by coated epoxy polyester resin and have a 30-mm high back lip to prevent accidents.



	Width		Depth	Height
Unit: mm	900 / 1200 / 1500 / 1800		750 / 900	900 / 750
Front View				
Model	<b>BE-900</b>	<b>BE-1200</b>	<b>BE-1500</b>	<b>BE-1800</b>



## LABORATORY BENCHES

### Island Laboratory Bench:

Central laboratory bench based on a modular of self-supporting steel structure made of cold rolled steel pipe, rectangular and circular, treated and covered with Epoxy Polyester coated.

In the middle of the benches, a structure supports creating a space for the guiding of the mechanical, electrical, plumbing, connections with the required outlets, easily accessible via removable covers. The bench frames are provided with leveling legs.

Service galleries contain shelves and service spines that can be configured depending on the client's needs. These storage shelves are anchored to the vertical profiles and easily height adjustable. They are protected by coated epoxy polyester resin and have a 30-mm high back lip to prevent accidents.



	Width		Depth	Height
Unit: mm	900 / 1200 / 1500 / 1800		1500 / 1650	900 / 750
Front View				
Model	<b>IBE-900</b>	<b>IBE-1200</b>	<b>IBE-1500</b>	<b>IBE-1800</b>



## LABORATORY BENCHES

### Movable Benches:

Movable benches are C frame benches mounted on 4 wheels, 2 of them are lockable. This design creates a flexible working environment for the end user.







	Width	Depth	Height
Unit: mm	1200 / 1500	750 / 900	900
Front View			
Model	<b>MBE-1200</b>		<b>MBE-1500</b>



## H-FRAME SYSTEM

The H-frame system was designed to provide a heavy-duty modular solution with under-bench registrable cover panels for easy maintenance & disinfection. Our H-frames are extremely stable and can load up to 300 kg/m<sup>2</sup>.



	Width		Depth	Height
Unit: mm	900 / 1200 / 1500 / 1800		750 / 900	900
Front View				
Model	<b>HBE-900</b>	<b>HBE-1200</b>	<b>HBE-1500</b>	<b>HBE-1800</b>

## Corner Panel

The corner panel for benches is used to connect the benches at the corner, made of laminated steel with epoxy coated. It has a registrable cover panel for easy maintenance.





	Width	Depth	Height
Unit: mm	750 / 900	750 / 900	750 / 900
Front View			
Model	<b>CP-A75D   CP-A90D</b>		



### ON PLINTH SYSTEM

The on-plinth system was designed to provide a solution with under-bench registrable cover panels for easy maintenance & disinfection. Underbench cabinets will be used to support worktops, without the need to use C or H-type frames. Each underbench unit is fitted with leveling legs & covered with a plinth.



	Knee space Width	Depth	Height
Unit: mm	830 / 1130	750 / 900	900 / 750
Front View			
Model	ZBE-1200	ZBE-1500	




## Balance Bench

Anti-vibration or balance bench has a self-supporting metal separate structure with elements of isolation and a vibration-free plate. Ideal for setting up analytical balances and other sensitive measuring equipment.

### General Characteristics:

- For weighing sensitive measurements with analytical balances
- Bench frame with worktop and vibration-free plate installed with 2 internal separate structures
- Not suited for laboratory operations requiring the extraction of dangerous chemicals



	Height	Width	Depth
Unit: mm	900	900 /1200	750 / 900
Front View			
Model	BB-900	BB-1200	









WORKTOPS



Arias collaborates with a number of worktop providers to ensure that the end user gets the best quality for their laboratory surfaces.

## Worktop Types:

	Ceramic	Epoxy	Phenolic Laminate (HPL)	Stainless Steel
<b>Material</b>	Ceramic	Resin with a hardener	Paper & resin based	Stainless Steel
<b>Corrosive Chemical Resistance</b>	High with no color change	High with noticeable color change	High with minor color change	Moderate
<b>Lab Application</b>	Moderate duty	Heavy duty	Moderate duty	Moderate to heavy duty
<b>Temperature Resistance</b>	High	High	Moderate to High	High
<b>Sample Environment</b>	Chemical, Industrial	Chemical, Industrial	Chemical, Industrial, Medical, & biological	Biology, food & beverages
<b>Scratch Resistance</b>	High	Weak to Moderate	High	Weak
<b>Images</b>				



SINK UNITS



## SINKS IN WALL OR ISLAND POSITION

Sink units are highly essential pieces of furniture in laboratory operations. Sinks can be connected to the various service systems available or be installed as individual units.

The Top sink is available in many materials: HPL, polypropylene, ceramic, stainless Steel, & epoxy.





The sink bowls are available in many materials like: Ceramic, Stainless Steel, Epoxy and Polypropylene.

	Ceramic	Epoxy	PP	Stainless Steel
<b>Material</b>	Ceramic	Resin with a hardener	Polypropylene	Stainless Steel
<b>Corrosive Chemical Resistance</b>	High with no color change	High with noticeable color change	High with no color change	Moderate
<b>Lab Application</b>	Moderate duty	Heavy duty	Moderate duty	Moderate duty
<b>Temperature Resistance</b>	High	High	Moderate	High
<b>Sample Environment</b>	Chemical, Industrial	Chemical, Industrial	Chemical, Industrial, Medical, & biological	Biology, food & beverages
<b>Scratch Resistance</b>	High	Weak to Moderate	Weak	Weak
<b>Images</b>				



### SINKS IN WALL OR ISLAND POSITION



	Width		Depth	Height
Unit: mm	600 / 900 / 1200 / 1500		750	900
Front View				
Model	<b>FRE-600</b>	<b>FRE-900</b>	<b>FRE-1200</b>	<b>FRE-1500</b>



STORAGE UNITS





### MATERIAL TYPES

All our storage units are made from steel epoxy-coated or melamine with an ergonomic design to ensure the safety, durability & functionality of the units.



Design Characteristics	Components
Load-bearing capacity per shelf (kg)	30
Construction	Body 1 part Adjustable shelves Door opening at 270° hinges 4 Paired polyamides, 2 with brake 4 heights adjustable feet
Combination possibilities	Drawers and doors
Handle	Zamak with in satin nickel
Drawers	Soft-Closing, fully extensible drawers



### STEEL STORAGE UNITS

Entirely built of cold-laminated galvanized steel, for different uses with high-quality accessories that can sustain harsh laboratory environments:

- Underbench cabinets
- Overhead cabinets
- Tall storage cabinets
- Open shelves units



### MELAMINE STORAGE UNITS

Made of melamine, for different uses with high-quality accessories that can sustain harsh laboratory environments:

- Underbench cabinets
- Overhead cabinets
- Tall storage cabinets



### CABINETS DESIGN

Arias provides a wide variety of storage cabinets:

- On-Wheel Cabinets
- Suspended Cabinets
- Sliding Cabinets
- On-Plinth Cabinets
- Corner Storage Cabinets
- Overhead Cabinets
- Tall Cabinets

They are all available in different sizes and designs to meet the needs & requirements of end users. All our cabinets comply with the following standard: "UNE EN 14727:2006.













## ON WHEELS

Under-bench cabinets on wheels are the best choice to have high flexibility for your benches to reconfigure them as you may wish in the future once you add new equipment. They can have 1 door, 1 door / 1 drawer or 4 drawers. Cabinets come with 4 wheels, 2 of them with breaks. Each cabinet includes an internal adjustable shelf.

### • Cabinet on wheels for H900

	Width			Depth	Height
Unit: mm	450 / 600 / 900			500	782
Front View					
Model	<b>CAB4-1PR</b>	<b>CAB4-1P1CR</b>	<b>CAB4-4CR</b>	<b>CAB9-2PR</b>	<b>CAB9-2P1CR</b>

### • Cabinet on wheels for H750

	Width	Depth	Height
Unit: mm	450	500	632
Front View			
Model	<b>CAB4-1PR-H75</b>	<b>CAB4-1P1CR-H75</b>	<b>CAB4-3CR-H75</b>



### SUSPENDED CABINETS

Underbench suspended cabinets are an innovative design for an easy cleaning and more free space. Cabinets can have 1 door, 1 door/1 drawer or 4 drawers. Each cabinet includes an internal adjustable shelf.






### SLIDING CABINETS

Sliding cabinets are the new design provided by Arias to support the evolving lab needs and requirements. It is a solution between suspended and on wheels as it allows end users to choose the position of the cabinet under the bench without having to use wheels. They can have 1 door, 1 door/1 drawer or 4 drawers. Each cabinet includes an internal adjustable shelf.








## • Suspended/Sliding Cabinet for H900

	Width			Depth	Height
Unit: mm	450 / 600 / 900			500	782
Front View					
Model	<b>CAB4-1PS</b>	<b>CAB4-1P1CS</b>	<b>CAB4-4CS</b>	<b>CAB9-2PS</b>	<b>CAB9-2P1CS</b>






## • Suspended/Sliding Cabinet for H750

	Width	Depth	Height
Unit: mm	450	500	632
Front View			
Model	<b>CAB4-1PS-H75</b>	<b>CAB4-1P1CS-H75</b>	<b>CAB4-3CS-H75</b>




## ON-PLINTH CABINETS

Cabinet On-plinth is the classic bench system that does not require table frames. This type of cabinets holds the worktop by itself. It can have 1 door, 1 door/1 drawer or 4 drawers. Each cabinet includes an internal adjustable shelf.

	Width			Depth	Height
Unit: mm	450 / 600 / 900			500	800
Front View					
Model	<b>CAB4-1PZ</b>	<b>CAB4-1P1CZ</b>	<b>CAB4-4CZ</b>	<b>CAB9-2PZ</b>	<b>CAB9-2P1CZ</b>

## CORNER STORAGE CABINET

A corner storage cabinet is a type of furniture that is designed to fit into the corner of a room, maximizing space and allowing for easy access to the contents inside. It is typically used to store items such as equipment, glassware or other supplies. The advantage of a corner storage cabinet is that it can help to save space in a room that might otherwise be unused. With adjustable shelves.

	Width	Depth	Height
Unit: mm	750 / 900	750 / 900	750 / 900
Front View			
Model	<b>CP-A75D   CP-A90D</b>		






## OVERHEAD CABINETS

Wall mounted cabinets are used to provide the end user with more storage space all around the lab walls to free more working areas. All have an internal adjustable shelf.

There are 3 types of cabinets available:

- With glass doors
- With steel doors
- With no doors



	Width		Depth	Height
Unit: mm	600 / 900 / 1200		370	700
Front View				
Model	<b>WCAB6-1P</b> <b>WCAB6-1PV</b>	<b>WCAB9-ENT</b> <b>WCAB12-ENT</b>	<b>WCAB9-2P</b> <b>WCAB12-2P</b>	<b>WCAB9-COR</b> <b>WCAB12-COR</b>







## TALL STORAGE CABINETS

Tall cabinets can be used for different storage purposes (glassware, balances and small devices...). The tall storage is available in different designs to adjust to the end user needs.

- Tall storage with 2 glass doors with hinges in upper part and 2 blinded doors in lower part
- Tall storage with 2 glass sliding doors in upper part and 2 blinded doors in lower part.
- Tall storage with 2 blinded doors
- Tall storage with 2 glass doors
- Tall storage with open shelves in upper part and 2 blinded doors in lower part

All shelves are height-adjustable.

	Width		Depth	Height
Unit: mm	900 / 1200		460	1950
Front View				
Model	<b>MA9-2P</b> <b>MA12-2P</b>	<b>MA9-2PV</b> <b>MA12-2PV</b>	<b>MA9-2P2V</b> <b>MA12-2P2V</b>	<b>MA9-2PENT</b> <b>MA12-2PENT</b>



### COLOURS


**Epoxy coated steel cabinets and metal components:**

For the steel cabinets and our benches metal components, we use 3 standard colours.



**Melamine particle board cabinets:**

The melamine cabinets have MDF laminate panels of 18 mm thickness with a wide choice of colours on the coated surface.





Other Storage Units designs:



### OPEN SHELF

Made of epoxy-coated steel, for different uses with adjustable shelves with a load capacity of 80 KG for each LM.

Open lab shelving is a type of storage system that is commonly used in laboratories, schools, and other research facilities. These shelving units are designed to hold and organize a variety of laboratory equipment, materials, and supplies. They are called "open" because they do not have doors or other enclosures, allowing easy access to the items stored on the shelves.

In addition to organizing lab equipment, open lab shelving can be used to store chemicals, reagents, and other lab materials that are used in experiments and research.



	Height	Width	Depth
Unit: mm	1800	900 / 1200	500
Front View			
Model	<b>ESTAL-900</b>	<b>ESTAL-1200</b>	



### TALL CHEMICAL CABINETS

Internal volume: approx. 408 litres.

#### Accessories

Possibility to provide the cabinet with activated charcoal filters and fan.

Constructed entirely from pickled and electrogalvanized steel, coated with powdered epoxy acid-resistance, then passage through a 200° C. thermic tunnel. Painting thickness of no less than 60 m (micron).

Completely sealed external monolithic structure. New construction using internal and exterior double cases.

Included electrogalvanized steel basins with watertight bottoms coated with epoxy anti-acid.

- Doors opening to 180°
- Key lock
- Acid-proof hinges
- "Danger" warning stickers on the doors



Tall Chemical Storage		Sizes		
Width (mm)	600mm	900mm	1200mm	
Depth (mm)	500mm	500mm	500mm	
Height (mm)	1980mm	1980mm	1980mm	



SERVICE GALLERIES



At Arias we have a wide range of service galleries that can be customized as per client's requirements in order to ensure a safe & flexible working environment. The ergonomic design of the service galleries allows easy access to connection points, permitting the passage of pipes inside it with outlets at the required position, avoiding any running over at the working area.

### Service Gallery Types:

- Wall Gallery
- Island Gallery
- Horizontal Service Spine
- Shelves Models







## WALL GALLERY

We have 2 types of wall gallery:





- **With Service Spine**

Gallery for wall benches with service spine system with ergonomic design that allows easy access to connection points, permitting the passage of pipes inside it with outlets at the required position, avoiding any running over the working area.




The service spine could be complemented by 1 or 2 raw shelves facing the user. Built according to Standard EN-14135, these shelves are easily height-adjustable without need to have them dismantled.

Shelves are available in steel or in glass.



	Width		Depth	Height
Unit: mm	900 / 1200 / 1500 / 1800		400 / 650	2120
Front View				
Model	<b>GCA-90C</b> <b>GCA-90</b>	<b>GCA-120C</b> <b>GCA-120</b>	<b>GCA-150C</b> <b>GCA-150</b>	<b>GCA-180C</b> <b>GCA-180</b>



	Width	Depth	Height
Unit: mm	1200 / 1500 / 1800	400 / 650	2120
Front View			
Model	<b>GCB-120C</b> <b>GCB-120</b>	<b>GCB-150C</b> <b>GCB-150</b>	<b>GCB-180C</b> <b>GCB-180</b>



### • Without Service Spine

Gallery for wall benches without service spine system used to have services and outlets top mounted. The position of outlets has to be defined previously by the client.

The service spine could be complemented by 1 or 2 raw shelves. Built according to Standard EN-14135, these shelves are easily height-adjustable without need to have them dismantled. Shelves are available in steel or in glass.







## ISLAND GALLERY

We have 2 types of island gallery:

- **With Service Spine**




Gallery for island benches with service spine system with ergonomic design that allows easy access to connection points, permitting the passage of pipes inside it with outlets at the required position, avoiding any running over the working area. The service spine could be complemented by 1 or 2 raw shelves on both sides. Built according to Standard EN-14135, these shelves are easily height-adjustable without need to have them dismantled. Shelves are available in steel or in glass.



	Width		Depth	Height
Unit: mm	900 / 1200 / 1500 / 1800		400 / 650	2120
Front View				
Model	<b>GCA-90C</b> <b>GCA-90</b>	<b>GCA-120C</b> <b>GCA-120</b>	<b>GCA-150C</b> <b>GCA-150</b>	<b>GCA-180C</b> <b>GCA-180</b>



## SERVICE GALLERIES

	Width	Depth	Height
Unit: mm	1200 / 1500 / 1800	400 / 650	2120
Front View			
Model	<b>GCB-120C</b> <b>GCB-120</b>	<b>GCB-150C</b> <b>GCB-150</b>	<b>GCB-180C</b> <b>GCB-180</b>



### • Without Service Spine

Gallery for island benches without service spine system used to have services and outlets top mounted. The position of outlets has to be defined previously by the client.

The service spine could be complemented by 1 or 2 raw shelves on both sides. Built according to Standard EN-14135, these shelves are easily height-adjustable without need to have them dismantled.




Shelves are available in steel or in glass.



### HORIZONTAL SERVICE SPINE

Laboratory trunking services supply at benches made of laminated steel and protected by epoxy paint. With registrable cover panels for easy maintenance. Modular system that can make different configurations. Integration of all service outlets including sockets and multiple connectors for information technology.



	Width	Depth	Height
Unit: mm	900 / 1200 / 1500	80	166
Front View			
Model	<b>CEM-90</b>	<b>CEM-120</b>	<b>CEM-150</b>



### SHELVES MODELS

All gallery systems could be complemented by 1 or 2 raw shelves for one or both sides. These are anchored to the vertical profiles with steel zinc coated injection brackets coated in epoxy-polyester resin with a metallic finish. These shelves are easily height-adjustable. They are protected by coated epoxy polyester resin; they also have a 30mm high back lip to prevent accidents.

Shelves are available in 2 types:


- **Steel**

Steel shelves protected by coated epoxy polyester resin.

- **Glass**

Tempered glass used on our shelves is 3 times stronger than regular glass. In addition to its strength, it also has excellent chemical resistance and light penetration to prevent and detect contaminations and leaks of reagent.






## WALL PARTITIONS & SERVICE COLUMN





## WALL PARTITIONS & SERVICE COLUMN



**Wall Partitions** are a growing demand in research laboratories. A movable table can be added to the wall partition and use its service spine providing a maximum level of flexibility and adaptability in the laboratory.

The wall service spine includes modular and interchangeable panels, that can be adapted to any need.

Wall partitions include back panels to have an easy access to the services and registrable cover panels for easy maintenance.

**Available in two options:**

- Full steel
- Full steel with glass



## WALL PARTITIONS & SERVICE COLUMN


A **Service column** has the main function of service connections from the ceiling and/or floor to the laboratory benches. Providing services to users, through gallery or horizontal electrical trunking systems, and may contain services in the panels, (gases, water, electricity and data).

The modular system's adaptability to different dimensions is possible through its flexibility and evolution, configurations and needs. It provides the opportunity to be central for central benches, mural for wall benches and finally they can have the total height of the room, with mural fixation and bench fixation. A service column is equipped with registrable cover panels, which facilitates easy maintenance and, if necessary, the future connections of services.

### Services supplied from the ceiling or from the floor for:

- Laboratory devices on mobile tables or under-bench positions
- Floor-mounted laboratory devices
- Version with one or two sides
- Design versions for genetical engineering areas
- Additional storage area through the connection of service modules with reagent repositories.



	Width	Depth	Height
Unit: mm	330 / 900 / 1200 / 1500	150	3000
Front View			
Model	<b>BAJG-330</b>	<b>TAB-90 TAB-120 TAB-150</b>	<b>TABV-90 TABV-120 TABV-150</b>



ACCESSORIES



## ACCESSORIES

At Arias, we are using a set of high quality accessories, that can withstand the harsh environment of the lab to ensure the product safety, durability and quality.

### LABORATORY TAPS

A large choice of laboratory taps that comply with the most demanding regulations for taps and fittings (DIN12918), for handles (EN13792) and for nozzle fittings (DIN12898).



### LABORATORY PRESSURE REGULATORS

For high purity gases (like carbon dioxide, hydrogen, helium & nitrogen), a large choice of regulators with gauge, inlet at rear, outlet downwards is provided.

A great variety of assembly possibilities in laboratory furniture with the required Gas type specific colour indication labels according to DIN13792.

May be combined with inlet shut-off valves, wall connectors, metering valves, different gauges and diverse accessories.

### ELECTRICAL SOCKETS

Electrical socket supports compatible with containers (European and British standard). They are available in different colors.






### PEG BOARDS

Made with High Density PP material, our pegboard standard size is L400\*W108\*H550mm. With droplet water rack, single faced with a drainage hole in the middle of one-sided bottom tray for water drainage. It has removable drip stick with locking function, easy to use.

**Detachable peg:** Convenient for use

**Installation:** Could be panel mounted or deck mounted

**Colour:** Available in grey and black



### SAFETY SHOWER

Made by galvanized steel piping with brass fittings and anti-corrosive coating in high visibility yellow or red colour. Recommended water inlet pressure 2 bar/Min. 1.5 bar/Max. 8 bar.

Shower head in ABS plastic and could include an eyewash as an additional option. It has two high flow aerated water sprays at low pressure with automatic opening anti-dust cover, and a push handle for operation.



### EYE WASH

Can be a double or 1 single mouth eyewash. The eye wash body is made by solid brass. Coating: High gloss epoxy powder, resistant to most chemicals, UV fading & heat. Spray head: soft rubber, softened water spray for eye protection. Dust Cup: PP, removed automatically when eye wash is operated. PVC pipe with stainless steel wire braided, wrapped by PE pipe for rust and leakage relief. Max. working pressure: 6 Bar.





### ARM EXHAUST

Local extractor arm hood. An ergonomic and space saving solution for simple and light applications in the Laboratory. Ducting and exhaust fan of this arm comes separate and are quoted separately.



### LABORATORY CHAIR

Ideally suited for tall counter height work areas.

#### Specifications:

- Variety of colours available.
- Chroming steel footrest ring.
- Nylon ESD 5 wheels.
- ESD Polyurethane foaming chair.

### LABORATORY STOOL

Ideally suited for desk work areas.

#### Specifications:

- Variety of colours available.
- Height adjustment by gas spring (42-54 cm)
- Artificial leather covered seat.
- Plastic coated skeleton.
- 5 pieces rotating castors.















## PRODUCT REFERENCE LIST



## PRODUCT REFERENCE LIST






### Sink Unit



	Width	Depth	Height	
Unit: mm	600 / 900 / 1200 / 1500	750	900	
Front View				
Model	FRE-600	FRE-900	FRE-1200	FRE-1500




### Movable Cabinet



	Width		Depth	Height	
Unit: mm	450 / 600 / 900		500	782	
Front View					
	Model	CAB4-1PR	CAB4-1P1CR	CAB4-4CR	CAB9-2PR

### Movable Cabinet for H750



	Width	Depth	Height
Unit: mm	450	500	632
Front View			
Model	CAB4-1PR-H75	CAB4-1P1CR-H75	CAB4-3CR-H75






### Suspended / Sliding Cabinet



g Cabinet

	Width	Depth	Height
Unit: mm	450 / 600 / 900	500	782

Front View



Model	CAB4-1PS	CAB4-1P1CS	CAB4-4CS	CAB9-2PS	CAB9-2P1CS
-------	----------	------------	----------	----------	------------



## PRODUCT REFERENCE LIST

### Suspended / Sliding Cabinet for H75

	Width	Depth	Height
Unit: mm	450	500	632



Front View



Model	CAB4-1PS-H75	CAB4-1P1CS-H75	CAB4-3CS-H75
-------	--------------	----------------	--------------

### On-Plinth Cabinet

	Width	Depth	Height
Unit: mm	450 / 600 / 900	500	880



Front View



Model	CAB4-1PZ	CAB4-1P1CZ	CAB4-4CZ	CAB9-2PZ	CAB9-2P1CZ
-------	----------	------------	----------	----------	------------

### Wall Bench

	Width	Depth	Height
Unit: mm	900 / 1200 / 1500 / 1800	750 / 900	900 / 750



Front View



Model	BE-900	BE-1200	BE-1500	BE-1800
-------	--------	---------	---------	---------

### Island Bench

	Width	Depth	Height
Unit: mm	900 / 1200 / 1500 / 1800	1500 / 1650	900 / 750



Front View



Model	IBE-900	IBE-1200	IBE-1500	IBE-1800
-------	---------	----------	----------	----------



## PRODUCT REFERENCE LIST

### Movable Bench

	Width	Depth	Height
Unit: mm	1200 / 1500	750 / 900	900



Model **MBE-1200** | **MBE-1500**

### H-Frame Bench

	Width	Depth	Height
Unit: mm	900 / 1200 / 1500 / 1800	750 / 900	900



Front View



Model **HBE-900** | **HBE-1200** | **HBE-1500** | **HBE-1800**

### On-Plinth Bench

	Knee space Width	Depth	Height
Unit: mm	830 / 1130	750 / 900	900 / 750



Front View



Model **ZBE-1200** | **ZBE-1500**

### Balance Bench

	Width	Depth	Height
Unit: mm	900 / 1200	750 / 900	900



Front View



Model **BB-900** | **BB-1200**








## PRODUCT REFERENCE LIST




### Gallery System for Wall and Central Tables



	Width	Depth	Height
Unit: mm	1200 / 1500 / 1800	400 / 650	2120
Front View			
Model	GCB-120C GCB-120	GCB-150C GCB-150	GCB-180C GCB-180




### Wall Partitions



	Width	Depth	Height
Unit: mm	330 / 900 / 1200 / 1500	150	3000
Front View			
Model	BAJG-330	TAB-90 TAB-120 TAB-150	TABV-90 TABV-120 TABV-150

### Service Spine



	Width	Depth	Height
Unit: mm	900 / 1200 / 1500	80	166
Front View			
Model	CEM-90	CEM-120	CEM-150



## PRODUCT REFERENCE LIST

### Corner Panel

Unit: mm

Width	Depth	Height
750 / 900	750 / 900	750 / 900



Front View



Model

CP-A75D  
CP-A90D

### Corner Storage

Unit: mm

Width	Depth	Height
750 / 900	750 / 900	750 / 900



Front View



Model

CP-A75D  
CP-A90D







## CERTIFICATIONS



## CERTIFICATIONS

### ARIAS LABS SYSTEM CERTIFICATIONS

Our company meets the highest European standards, quality and safety are our absolute priority.



#### TEST CERTIFICATE n° 230.Z.1608.479.EN.01

References: 1606142-01 – 1707049-01 – Cj

**PRODUCT:** Laboratory module hung unit,  
Code: CAB4C-1PD1C.  
Additional code module: CAB4C-1PD,  
CAB4C-1PI, CAB4C-1PI1C, CAB4C-4C



**COMPANY:** ARIAS LABS SYSTEM, S.L  
St/ FUENTE ALAMO, 16  
28400 COLLADO VILLALBA (MADRID)  
Phone.: +34 808 51 55 24  
VAT: B 87540787

**TEST:** Compliance with the following standards:  
"UNE EN 14727:2006"  
Laboratory furniture - Storage units for laboratories - Requirements and test methods.

**RESULT:** The model tested satisfactorily fulfills the especifications for the standards applied in the following sections:

TESTS	RESULT
- Sect. 6 - General safety requirements	CORRECT
- Sect. 7.1.1 - Stable position of the shelves	CORRECT
- Sect. 7.1.2 - Deflection of shelves	CORRECT
- Sect. 7.2 - Strength test of shelf supports	CORRECT
- Sect. 7.3.2 - Strength of pivoted doors	CORRECT
- Sect. 7.3.4 - Durability of pivoted doors with hinges (n= 50.000 cycles)	CORRECT
- Sect. 7.5.2 - Slam open of drawers	CORRECT
- Sect. 7.5.3 - Strength of the drawers	CORRECT
- Sect. 7.5.4 - Durability of the drawers (n= 50.000 cycles)	CORRECT
- Sect. 7.5.5 - Displacement drawer bottom	CORRECT
- Sect. 7.7.2 - Static load on top surfaces of storage units	CORRECT

Paterna, September 14, 2017

**AIDIMME®**  
Signed: José Emilio Nuévalos  
Head of Furniture Laboratory

This certificate only refers to the samples tested by the AIDIMME laboratory.  
The particular results of the tests are described in technical report N° 230.I.1709.522.ES.01 dated on 07/24/2017

AIDIMME is a member of INNOVAWOOD, The European Network of Research and Training for the Forest, Wood and Furniture Industry, among whose members are: BRE-CTTC (United Kingdom), COSMOB (Italy), DTI (Denmark), FCBA (France), ITD (Poland), SHR (Holland), SP Trätek (Sweden), TRADA-FIRA (United Kingdom), University of Zagreb (Croatia), WKI (Germany).

AIDIMME: METAL-PROCESSING, FURNITURE, WOOD AND PACKAGING TECHNOLOGY INSTITUTE  
Tax address: Parque Tecnológico - Calle Benjamín Franklin, 13  
CIF: E5G46261590 - 46980 PATERNA (VALENCIA) ESPAÑA  
Phone: +34 96 136 60 70 - Fax: 96 136 61 89  
Registered office: Parque Tecnológico - Avda. Leonardo Da Vinci, 36  
46980 PATERNA (VALENCIA) ESPAÑA  
Phone: +34 96 131 85 97 - Fax: 96 001 54 46

aria@aidimmes.com  
www.aidimmes.com



## CERTIFICATIONS

### ARIAS LABS SYSTEM CERTIFICATIONS

ASHRAE110-1995:



ANSI/ASHRAE 110-1995Flow Containment & Instrumentation

**CERTIFICATE OF TYPE TESTING IN  
ACCORDANCE WITH ASHRAE 110-1995**

**CERTIFICATE & REPORT NO: INV/ASHRAE110/872**  
**ISSUE DATE: 25th September 2018**

**Fume Hood Supplier:**  
**Arias Labs System**  
**Hermanos Lumiere 12**  
**28500 Arganda**  
**Madrid**  
**SPAIN**



**Fume Hood Type:**  
**1.5m wide bench-type vertical-sash**  
**External Dimensions:**  
*Height = 2400mm*  
*Width = 1500mm*  
*Depth = 900mm*  
**Internal Dimensions:**  
*Height = 1060mm*  
*Width = 1400mm*  
*Depth = 660mm (wall-sash)*  
*600mm (baffle-sash)*  
**Baffle Gap Dimensions:**  
*Top gap: 25mm*  
*Bottom gap: 200mm*

**Test Opening:**  
**Width: 1200mm**  
**Height: 500mm**  
  
**Fume Hood Flow:**  
**Flow rate: 0.340m<sup>3</sup>/s**  
**Face velocity: 0.52m/s**  
**Press drop: 140Pa (at 250mmD)**  
  
**Fume Hood Containment:**  
**Position P1: <0.020ppm**  
**Position P2: <0.020ppm**  
**Position P3: <0.020ppm**  
**Edge scan: <0.020ppm**

This is to certify that the fume hood described above has been type-tested in accordance with ASHRAE 110-1995 and resulted in performance characteristics given in the corresponding test report.

**Tested and Certified by: Dr A F Bicen**

**Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR**  
**T: +44 (0) 1582 461666 F: +44 (0) 1582 460077 E: info@invent-uk.com**






## CERTIFICATIONS

### ARIAS LABS SYSTEM CERTIFICATIONS


EN14175:



EN 14175Flow Containment & Instrumentation

**CERTIFICATE OF TYPE TESTING IN  
ACCORDANCE WITH EN 14175, PART 3**  
**CERTIFICATE & REPORT NO: INV/EN14175/871**  
**ISSUE DATE: 25th September 2018**

**Fume Cupboard Supplier:**  
**Arias Labs System**  
**Hermanos Lumiere 12**  
**28500 Arganda**  
**Madrid**  
**SPAIN**




**Fume Cupboard Type:**  
**1.5m wide bench-type vertical-sash**  
*External Dimensions:*  
*Height = 2400mm*  
*Width = 1500mm*  
*Depth = 900mm*  
*Internal Dimensions:*  
*Height = 1060mm*  
*Width = 1400mm*  
*Depth = 660mm (wall-sash)*  
*600mm (baffle-sash)*  
*Baffle Gap Dimensions:*  
*Top gap: 25mm*  
*Bottom gap: 200mm*

**Test Opening:**  
**Width: 1200mm**  
**Height: 500mm**  
**Fume Cupboard Flow:**  
**Flow rate: 0.340m<sup>3</sup>/s**  
**Face velocity: 0.518m/s**  
**Press drop: 140Pa (at 250mmD)**  
**Fume Cupboard Containment:**  
**Inner-plane containment:**  
**C1: <0.010ppm**  
**Outer-plane containment:**  
**C2, C3, C4, C5: <0.010ppm**  
**Robustness of containment:**  
**CR: 0.011ppm**

This is to certify that the fume cupboard described above has been type tested in accordance with Part 3 of EN 14175, in compliance with the requirements of Part 2 and with reference to Part 1, and resulted in performance characteristics given in the corresponding test report.

**Tested and Certified by: Dr A F Bicen**



**Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR**  
**T: +44 (0) 1582 461666 F: +44 (0) 1582 460077 E: info@invent-uk.com**



## CERTIFICATIONS

### ARIAS LABS SYSTEM CERTIFICATIONS



**TEST CERTIFICATE n° 230.Y.1605.282.EN.01**

Reference: 1604013-01 – C-1



**PRODUCT:** Table for laboratory model:  
"ARIAS / BF9-12-D / BASIC FRAME 1200"

**COMPANY:** ARIAS LABS SYSTEM, S.L  
St./ FUENTE ALAMO, 16  
28400 COLLADO VILLALBA (MADRID)  
Phone: 608 51 55 24  
VAT: B 87540787

**TEST:** "UNE EN 13150:2001"  
WORKBENCHES FOR LABORATORIES. DIMENSIONS, SAFETY  
REQUIREMENTS AND TEST METHODS.

**RESULT:** The model tested satisfactorily fulfills the specifications for the standard applied in  
the following tests carried out:

TEST	RESULT
- Sect. 4 – Recommended Dimensions	CORRECT
- Sect. 5 – General safety requirements	CORRECT
- Reflection of the surface according to (ISO 2813) to 60°	13,6
- Sect. A.3.1 - Horizontal Static Load Test (Fh= 600 N)	CORRECT
- Sect. A.3.2 - Vertical Static Load test (Fv= 2000 N)	CORRECT
- Sect. A.3.5 - Stability Test Versus Horizontal Impact (M=50 kg, h = 40 mm)	STABLE
- Sect. A.3.6 - Stability under vertical load (Fv= 1000 N)	STABLE

Paterna, 29<sup>th</sup> April 2016



Signed: José Emilio Nuévalos  
Head of Furniture Laboratory

This certificate only refers to the samples tested by the AIDIMA laboratory.

The particular results of the tests are described in technical report ref.: 1604013-01 dated on 29/04/2016.

AIDIMA is a member of INNOVAWOOD, The European Network of Research and Training for the Forest, Wood and Furniture Industry, among whose members are: CATAS SPA (Italy), COSMOB (Italy), CTIB-TCHN (Belgium), DTI (Denmark), ELKEDE (Greece), EMPA (Switzerland), FCBA (France) HFA (Austria), IHD (Germany), ITD (Poland), SHR (Holland), SP Trätekt (Sweden), TRADA-FIRA (United Kingdom), University of Zagreb (Croatia), VTT (Finland), WKI (Germany).

Parque Tecnológico - Calle Benjamín Franklin, 13  
Avenida nº 50 - 46100 PATERNA (Valencia) ESPAÑA  
Tel: 96 136 60 70 - Fax: 96 136 61 85E-mail: aidima@aidima.es

[www.aidima.es](http://www.aidima.es)







Calle N3, Parcela C  
Polígono industrial el Bobalar,  
46130 Massamagrell, Valencia  
Tel: +34 960 013 119 | +34 634 652 807

**[www.ariaslabssystem.com](http://www.ariaslabssystem.com)**